

# EXHIBIT 3

# NasalGuard®

AIRBORNE PARTICLE BLOCKER

## Product Facts

NasalGuard® Airborne Particle Blocker is a clear gel, applied outside the nose, that uses electrostatic technology to reduce the inhalation of airborne particles before entering your nasal passages.

## Uses

Appropriate for everyday use or prior to exposure to harmful particles including contaminated air, virus-sized particles, mold, pollen, dust mites, and pet dander. Safe for everyone including pregnant women, nursing mothers, and children.

## Warning

Avoid contact with eyes. Rinse with water if accidental eye contact or skin irritation occurs. Discontinue use if any skin irritation or rash occurs.

## Directions

For external use only. Apply before exposure to airborne contaminants.



1 Begin by squeezing 1 to 2 pin-sized drops on your pointer finger.

2 Spread the gel between your pointer finger and thumb to ensure an even application.

3 Apply directly around the nostrils and above the upper lip.

Apply every 4 to 6 hours or as often as needed. Reapply if nose or face becomes wet.

READ AND SAVE THIS CARD FOR REFERENCE

## Questions or Comments?

Visit [www.NasalGuard.com](http://www.NasalGuard.com) or call 1-855-NASALGUARD (1-855-627-2548)

Distributed by:  
NasalGuard LLC.  
Somerville, NJ 08876, USA

US Patent # US 6,844,005 B2 US 8,163,802 B2  
Other Patents Pending

© NasalGuard LLC 2017

LOT NsGd01

Made in USA

